
# Ultrasonic Cleaner

Enzyme based concentrate for the pre-cleaning and decontamination of the non-porous surfaces of hand instruments and forceps, rotary instruments, mirrors etc. Designed for use in the Ultrasonic bath. Specially formulated for use in the Medical and Dental environment.

**Directions for Use:** Add 30ml ± 1ml of the concentrate to 970ml of tap water in the ultrasonic bath. (3% solution). Immerse the instruments ensuring that all surfaces are covered by the solution, add hollow objects in such a way that air escapes. Run the ultrasonic cleaner for a minimum of 15 minutes and according to the manufacturers instructions. After immersion, remove the instruments and rinse them thoroughly prior to terminal sterilisation or high level disinfection. Prepare a fresh solution daily or sooner if solution is very dirty. For professional use only.

**Composition:** Aqua, Propan-2-ol, Didecylidimethylammonium chloride, Sodium n-lauroylsarcosinate, Subtilisin, Dye.

**Warning.** H226: Flammable liquid and vapour. H315: Causes skin irritation. H319: Causes serious eye irritation. H336: May cause drowsiness or dizziness. P302+350: IF ON SKIN: Gently wash with plenty of soap and water. P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P280: Wear protective gloves/protective clothing/eye protection/face protection. P261: Avoid breathing dust/fumes/gas/mist/vapours/spray. P304+340: IF INHALED: Remove person to fresh air and keep comfortable for breathing. **Storage.** P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P403+P233: Store in a well-ventilated place. Keep container tightly closed. P235: Keep cool. **Contraindications.** Incompatible materials: Aluminium, copper. Not recommended for plastics and rubbers. **Disposal.** Do not discharge concentrate into drains or rivers.

UK	Warning.
BG	ВНИМАНИЕ.
CZ	Varování.
DE	Achtung.
DK	Warnung.
EE	Hoiatus.
ES	Atención.
FI	Varoitus.
FR	Attention.
GR	Προσοχή.
HR	Upozorenje.
HU	Figyelem.
IT	Attenzione.
LT	Ispėjimas.
LV	Bīdinājums.
NL	Waarschuwing.
NO	Advarsel.
PL	Uwaga.
PT	Atenção.
RO	Atenție.
SE	Varning.
SI	Pozor.
SK	Pozor.

Topdental (Products) Ltd.  
12 Ryefield Way, Silsden,  
W. Yorks, BD20 0EF,  
England, UK.  
Tel +44 (0)1535 652 750  
www.topdental.co.uk

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