


## Instrument Disinfectant Powder with Colour Change


Date of revision: 20/11/24

Available in 700g


**Intended Use / Intended Purpose:** Disinfectant powder for use in the reprocessing of all reusable surgical, medical and dental instruments prior to sterilisation. For use in the instrument soaking bath.




1. Dissolve 1 level measure of powder (20g ± 2g) in 980ml of warm tap water (30 - 40°C).




2. Place both water and powder into the instrument bath




3. Mix for 60 seconds. The solution will turn red.



4. Completely immerse the instruments into the solution ensuring that no air is left within hollow instruments.



4. Leave for 30 minutes.



5. After disinfection, rinse instruments thoroughly with RO or distilled water prior to further processing and sterilisation. The solution will remain effective for 8 hours. Change to fresh solution daily or earlier should physical dirt in the solution become evident. We recommend preparing fresh solution at the start of the working day. When the solution is colourless it is no longer active. See product labels and MSDS for further information.

**Indications / Efficacy (30 mins):** Bactericidal including MRSA. Yeastcidal. Sporicidal. Viruses to include: Non-enveloped viruses: Murine norovirus, Adenovirus 5. Enveloped viruses: Hepatitis B virus (HBV), Hepatitis C virus (HCV), Hepatitis Delta virus (HDV), HIV, Influenza virus, Coronavirus, Herpesviridae, Filoviridae (Ebola), Rubella, Rabies, Measles, Poxviridae.

**Tested to Standards:** EN17126, EN 17111, EN14885, EN14562, EN14561, EN 14476, EN13727, EN13624.

**Composition:** TAED, Sodium Percarbonate, Citric Acid Monohydrate, Sodium Alkyl Aryl Sulphonate, Sodium Sulphate Anhydrous, Tetrasodium (1-hydroxyethylidene) bisphosphonate, Sodium Nitrite, Dye.

**Contraindications.** Aluminium, Copper.

**Disposal.** Do not discharge into drains or rivers.

**Intended Users.** For professional use only. The products are used within the dental surgery by the dentist or dental assistant as part of their general disinfecting routines.

**Patient Population:** Dental surgery patients, both adults and children.

Any serious incident related to this medical device must be reported to the manufacturer and the competent authorities.

**Danger.** H318: Causes serious eye damage. H272: May intensify fire; oxidiser. H411: Toxic to aquatic life with long lasting effects. P280: Wear eye protection, protective gloves. P305+P351+P338+P310: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a doctor. P391: Collect spillage. **Storage.** P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P220: Keep away from clothing and other combustible materials. Store in a well-ventilated place. Keep container tightly closed. Keep dry.



**Danger.**



For additional languages see [www.topdental.co.uk](http://www.topdental.co.uk)



Manufactured in the UK by:  
Topdental (Products) Ltd.  
Unit 1, Holmfield Industrial Estate,  
Halifax, W. Yorks, HX2 9TN. UK  
Tel +44 (0)1535 652750  
[www.topdental.co.uk](http://www.topdental.co.uk)



DIRECTA AB P.O. Box 723,  
Finvids väg 8 SE-194 27  
Upplands Väsby, Sweden



UFI: QJ10-F0AF-F002-H9YP