



ECO-XL PLASTIC FREE WIPES 240 x 200mm

- Plastic free wipe
- Plant based fibres
- Renewable resources
- Reduced carbon footprint
- Alcohol free

V-WIPE ZERO Surface Disinfectant Wipe BIODEGRADABLE[†]



2 in 1
DISINFECTS
& CLEANS



KEY FEATURES

V-WIPE Zero ECO-XL wipes are plastic free wipes made from renewable and sustainable natural viscose fibres. A modern and eco friendly alternative to traditional wipes which are made from Polypropylene or Polyester. Developed for use in clinical environments, on non invasive medical devices and general high traffic non-porous surfaces.

The new, exclusive wipe substrate material is constructed from plant based viscose fibres which are 100% renewable and sustainably sourced. Wipes are biodegradable[†] and totally plastic free.

The new flowpack style packaging helps to reduce the amount of plastic consumed by removing the dispenser tub. Each pack includes an easy open and sturdy stay fresh lid. This ensures wipes remain moist and ready for the next application.

The alcohol free formulation is ideal to use on the majority of wipeable surfaces even those that are sensitive to alcohol, including dental chairs*

WHAT YOU NEED TO KNOW

- Biodegradable[†]. Plastic free wipes
- 100 extra large wipes
- Contact time from only 60 seconds
- Dual registration; Medical Device and Biocide
- Phase II Step I & II tested
- Safe to use on sensitive surfaces*
- 2 in 1; Cleans and disinfects
- Subtle apple aroma

Bossklein ECO-XL contains a synergistic blend of two powerful biocides, is free from PHMB and contains detergent derived from renewable resources.

The new formula takes advantage of modern, safer and greener chemical components to create a more simple and effective product for the end user.

V-WIPE Zero ECO-XL wipes are ideal to use where time is critical, such as in between patients and in high traffic areas.

COMPOSITION

Contains Aqua, Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/BKC (C12-16)) (68424-85-1) 0.48% w/w (4.8g/kg), Didecylmethylammonium Chloride (7173-51-5) 0.22% w/w (2.2g/kg), Perfumes.

INSTRUCTIONS FOR USE

Using sufficient wipes for each application. Wipe the clean surfaces to be disinfected so that they are completely moistened and let the disinfectant take action over 60 seconds (to achieve full efficacy leave for 2 minutes). Allow surface moisture to dry or wipe with clean tissue. Contains detergent. For professional use only.

Use biocides safely. Always read the label and product information before use.

[†] Wipes are made from biodegradable plant fibres

*Always test an inconspicuous area first See website for full details.
Do not use on copper or aluminium.



Complete hygiene solutions

DirectaDentalGroup
DIRECTA | Ispol dental | **osana parkell** | **KOHLER RONWIG** | **polydentic**
TrollDental | Sendoline | **physics** | **CONTACTE** | **ceramix**

Topdental (Products) Ltd
Unit 1, Holmfield Industrial Estate
Halifax, West Yorkshire HX2 9TN
Tel: +44 (0)1535 652 750 www.topdental.co.uk

Surface Disinfectant Wipe BIODEGRADABLE[†]

ECO-XL

- Plastic free wipe material
- Made from renewable resources
- Extra large wipe 240 x 200mm
- Alcohol free formulation
- Pack of 100 wipes
- Safe to use on most surfaces**

BIODEGRADABLE[†] MATERIAL

Unlike most competitor wipes which are made from polypropylene, Bossklein ECO-XL wipes are constructed from 100% natural plant based fibres. After use the plastic free, biodegradable[†] wipes will naturally breakdown overtime.

The construction and large size of the 40gsm, plant based, viscose wipe material is extra durable and ideal for clinical disinfection. The plastic free substrate is designed to help transfer the alcohol free disinfectant solution on to the surface in question.



PLANT
BASED WIPE



ECO-XL
PLASTIC FREE WIPES
240 x 200mm

MICROORGANISM	CONTACT TIME
Bactericidal	
BS EN 16615:2015	
BS EN 13727:2012 + A2:2015	
Pseudomonas aeruginosa	2 min
Staphylococcus aureus	2 min
Enterococcus hirae	2 min
Yeasticidal	
BS EN 16615:2015	
BS EN 13624:2021	
Candida albicans	60 sec
Enveloped viruses	
BS EN 14476:2013 + A2:2019	
BS EN 16777: 2018	
Vaccinia virus	60 sec
Hepatitis B virus (HBV)	60 sec
Hepatitis C virus (HCV)	60 sec
Hepatitis Delta virus (HDV)	60 sec
HIV	60 sec
Influenza virus	60 sec
Coronavirus	60 sec
Herpesviridae	60 sec
Filoviridae (Ebola)	60 sec
Rubella virus	60 sec
Rabies virus	60 sec
Measles virus	60 sec

MICROBIOLOGICAL EFFICACY

Bactericidal, Yeasticidal and Virucidal against all enveloped viruses including Coronavirus, Hepatitis B and C and HIV. This highly effective formula has been rigorously tested to real world EN standards for proven effectiveness including Phase II Step I and II.

The new alcohol free Bossklein Zero surface disinfectant is effective against all listed (opposite) micro-organisms. This ensures complete peace of mind when it comes to clinical surface disinfection.

The following test standards have been achieved: EN16777, EN16615, EN14885, EN14476, EN13727, EN13624.

VARIANTS AND PACKAGING

BOS8001 Flowpack of 100 ECO-XL Wipes x 9 per box



BACTERICIDAL



YEASTICIDAL



VIRUCIDAL*



[†] Wipes are made from biodegradable plant fibres

* Virucidal against all enveloped viruses

** Always test an inconspicuous area first. Do not use on copper or aluminium

