Bactericidal, Yeasticidal, All Enveloped Viruses inc Vaccinia, Hepatitis, HIV, Measles, Corona, +more Adenovirus 5, Norovirus,

2%-20<sup>ml</sup>/<sub>1L</sub>

 $15_{min}$ 



**DAILY**asp

Aspirator Cleaner & Disinfectant

Concentrated 2% (20ml per L) Low foaming non abrasive 15 min contact time

For everyday use with all types of dental aspirator / suction systems.

**SD**activ Surface Disinfectant Spray

Bactericidal, Yeasticidal, All Enveloped Viruses, Adenovirus 5, Norovirus, Mycobactericidal, Tuberculocidal



V-WIPE Surface Disinfectant Wipe

Contains alcohol (63%) 60 second contact time Refill and reuse tubs Ready to use sprays and wipes I emon aroma



**ULTRA**bio Ultrasonic Instrument Cleaner

Concentrated 3% (30ml per L) Enzymatic 15 min contact time

For use inside an ultrasonic bath to clean manual instruments inc scalpels, forceps, luxators, etc.

3%-30<sup>ml</sup>/1L



**IMP**activ | **P** Impression Disinfectant

Concentrated 2% (20g per L) - 1 scoop 60 sec contact time

For use with all types of dental impression materials including orthodontic impressions.

Bactericidal, Yeasticidal, All Enveloped Viruses inc Vaccinia,

2%-20<sup>9</sup>/1L | 60<sub>sec</sub>



Dactiv **Instrument Disinfectant** 

Concentrated 2% (20ml per L) Plant based ingredient 15 min contact time

For use inside a soaking or ultrasonic bath to clean manual instruments inc scalpels, forceps, luxators, etc.

and Cleaner

Bactericidal, Yeasticidal, All Enveloped Viruses in 15 mins using 2% Adenovirus 5 & Norovirus in 30 min using 4%, Mycobactericidal & Tuberculocidal n 60 min using 4%

2%-20<sup>ml</sup>/1L



DactivQUANTUM

**Sporicidal Instrument Disinfectant and Cleaner** 

Concentrated 2% (20g per L) - 1 scoop

30 min contact time For use inside a soaking bath to clean manual

Bactericidal, Yeasticidal, All Enveloped Viruses Adenovirus 5 & Norovirus, Sporicidal

instruments inc scalpels, forceps, luxators, etc.

2%-20%1L 30min

CLASSIC inspection lamps, hand controls, etc

For use on clinical surfaces resistant to alcohol inc worktops, READY TO USE WIPES & SPRAY **60**<sub>sec</sub> **UNI**oil Universal Handpiece Oil Aerosol lubricant spray High penetrative ability Contains peppermint oil Contains silicone and anti-rust inhibitors For use on all types of dental handpiece after each use for effective lubrication.

**SD**activ<sup>ZERO</sup> Surface Disinfectant Spray

Bactericidal, Yeasticidal, All Enveloped Viruses, Adenovirus 5, Norovirus, Mycobactericidal, Tuberculocidal



Surface Disinfectant Wipe

Alcohol free 60 second contact time Refill and reuse tubs Ready to use sprays and wipes



**DRY** wipes **Highly Absorbent Dry Wipe** Strong even when wet

For use with surface disinfectant sprays

**60**sec Contains detergent for added cleaning power.

For use in clinical areas on sensitive surfaces inc the dental chair\*