Bossklein



NEW VS OLD

The recent pandemic highlighted not only the requirement for effective hand disinfection but the differences between high grade, professional solutions and everyday, off the shelf products.

New B70gel from Bossklein is a high performance, medical grade hand sanitizer that has been designed to replace the outgoing Bossklein Hand Disinfectant Gel.

With its high alcohol content, extensive effective microbiological spectrum, and the addition of skin moisturising properties provided by aloe vera - which leaves hands super soft even when used frequently. B70gel is the ultimate in usable defence against contact transmission of pathogens between people and surfaces.

Below we will compare the main headline differences between the old and the new Bossklein hand disinfection products.



HIGHLIGHTS

- Bossklein B70gel contains a higher alcohol content (70%)
- Full dermatological testing has been undertaken for B70gel
- Contains aloe vera for softer feeling hands even after frequent use, unlike other gels
- Now with the powerful antioxidant vitamin E to help nourish and protect the skin

DISINFECT **PROTECT** RFFRFSH







NEW

Hand Disinfectant Gel Name

Variants BOSCLG050: 50ml Tottle with clip

BOSCLG500: 500ml Bottle with pump

BOSCLG005: 5 Litre Container

Ethanol (64-17-5), Aqua, Propylene Glycol, Car-Composition

bomer, Polysorbate-20, Aminomethyl Propanol.

Alcohol content Ethanol 62% w/w (620g/kg).

Contact time 30 seconds

Alcoholic Aroma

Colour Clear - Dve free

Dermatologically tested No Contains Vitamin E No

Contains aloe vera

B70gel Hand Disinfectant Gel

BOS7007: 50ml Tottle with clip BOS7006: 500ml Bottle with pump BOS7005: 5 Litre Container with pump

Ethanol (64-17-5), Aqua, Glycerine, Carbomer, Aminomethyl propanol, Aloe Barbadensis Leaf Juice, Vitamin E Acetate,

Ethanol 70% w/w 70g/100g

30 seconds

Aloe vera + alcohol

Clear - Dve free

Yes Yes





Bossklein[®]

B7000 ALOE VERA

Hand Disinfectant Gel







NEW

30 sec 30 sec

OLD

Efficacy Spectrum	Bactericidal	30 sec	;
	Pseudomonas aeruginosa	30 sec	
	Staphyloccocus aureus	30 sec	
	Enterococcus hirae	30 sec	
	E. coli	30 sec	
	Yeasticidal	30 sec	;
	Candida Albicans	30 sec	
	Enveloped Viruses	30 sec	;
	Vaccinia Virus	30 sec	
	Hepatitis B virus (HBV)	30 sec	
	Hepatitis C virus (HCV)	30 sec	
	Hepatitis Delta virus (HDV)	30 sec	
	HIV	30 sec	
	Measles Virus	30 sec	
	Influenza Virus	30 sec	
	Coronavirus	30 sec	
	Herpesviridae	30 sec	
	Filoviridae (Ebola)	30 sec	
	Rubella Virus	30 sec	
	Rabies Virus	30 sec	

Test Standards	EN 14476, EN 13727, EN 13624,
	EN 1500

Poxviridae

Bactericidal
Pseudomonas aeruginosa
Staphyloccocus aureus

Staphyloccocus aureus	30 sec
Enterococcus hirae	30 sec
E. coli	30 sec
Yeasticidal	30 sec
Candida Albicans	30 sec
Enveloped Viruses	30 sec
Vaccinia Virus	30 sec
Hepatitis B virus (HBV)	30 sec
Hepatitis C virus (HCV)	30 sec
Hepatitis Delta virus (HDV)	30 sec
HIV	30 sec
Measles Virus	30 sec
Influenza Virus	30 sec
Coronavirus	30 sec
Herpesviridae	30 sec
Filoviridae (Ebola)	30 sec
Rubella Virus	30 sec
Rabies Virus	30 sec
Poxviridae	30 sec

Norovirus EN14776, EN 13727, EN 13624,

Non-Enveloped Viruses

Adenovirus 5

EN 1500

HIGHLIGHTS

B70gel is effective against the hard to kill, non enveloped viruses Adenovirus and Norovirus

30 sec

- Bossklein B70gel keeps the rapid contact time of only 30 seconds
- All enveloped viruses are killed in only 30 seconds, including Coronavirus
- Aloe vera aroma leaves hands smelling fresh and feeling great

DISINFECT
PROTECT
REFRESH

30 sec

30 sec



